

BASIC

Certification Guide


by Zengar






1






2



3

The Zengar Basic Certification Guide




THE BRAIN COMMUNICATES BY DETECTION OF CHANGE

The **Language** our body speaks is **Dynamical**, when we get an infection the CNS immediately **'detects the change'** this tells our body to do something about the change it noticed.

The brain's ability to create pathways and **adapt to a new set of change is essential for learning and 'surviving'** ...this is what makes it a **'Self-Organizing Dynamical System'**

4

THE BRAIN AND CENTRAL NERVOUS SYSTEM (CNS)
THE BRAIN IS INFINITELY CYCLING AND SELF-REGULATING




Detects Changes by using the DIFS

- Duration
- Intensity
- Frequency
- Shift

5


THE BRAIN IS GEARED TOWARDS THE GENERAL



- The brain is **inherently geared towards general (but not absolute)** self-regulation based on information
- Ex. No need to teach a child the 'Best' way to walk – through experimentation the brain discovers and develops internal balance and a natural gait unique to that individual...
- **Brain** processes and assesses feedback from its own muscles & muscular-skeletal system about its position and stability in space

6

The Zengar Basic Certification Guide



PROCESS VARIABILITY

- Infinitely cycling & self regulating
- No **ABSOLUTE** values that are inherently 'best' or the specific "amount" of activity required for specific functions.
- It always varies from **Brain to Brain** and it always varies with time for that same **Brain**

7




NEUROFEEDBACK IS...

Neurofeedback is sometimes referred to as EEG (Electroencephalogram) Biofeedback or Neurotherapy.

It is the process of **training the brain/CNS** (central nervous system) to learn to modify and control electrical activity through auditory and visual feedback.

8



NEUROPTIMAL® DYNAMICAL NEUROFEEDBACK®


- It is safe
- No need for a diagnosis
- It is training not treatment
- It is a set of proprietary Mathematical Transforms
- JTFA and NDS (Nonlinear Dynamical Systems®)
- Unique Time-Frequency Envelopes (TFE) or Targets

9

The Zengar Basic Certification Guide



10



Neurofeedback – 2 Basic Paradigms

Linear Non-Dynamical NF- process uses a software program that measures electrical activity in 2-ways

- Amplitude
- Frequency

• **Nonlinear Dynamical NF-** Uses a software program that utilizes 4 ways to detect the changes in electrical activity-

- Duration-
- Intensity-
- Frequency
- Shifts

11



LINEAR MODELS

- The **linear model** believes: they need to affect the strength (amplitude) of the signal (by **augmenting and inhibiting** a specific frequency).
- The frequency was *chosen by looking at pre-recorded activity under different conditions (eyes open, closed, reading, counting by 7's starting at 783) this is called a **QEEG**. (Quantitative Electroencephalogram).
- This **Linear approach** uses **thresholds** that are static bars that 'stop' the sound, thereby initiating the feedback process

12



DIAGNOSTICALLY AGNOSTIC


- 'Regular Training' Session- never changing for any reason
- It never matters what the diagnosis is the **NeuroOptimal®** system will adapt to any 'brain' no matter if they medicated, addicted to a substance, or brain damaged in any way.

13



THE NEUROPTIMAL® SYSTEM IS A MATHEMATICAL DANCE BETWEEN THE CLIENT AND THE NO SYSTEM

14




THE NEUROPTIMAL® SOFTWARE IS...

Constantly **Adapting** to the individual's electrical activity **256 times per sec**

Is specific to the **Individual on the NO System** at that moment in time **Changing** as the **Individual Changes (Dancing)**

15

The Zengar **Basic Certification Guide**



16





NEUROPTIMAL® DETECTS “SHIFTS IN EEG”

Shifts in EEG are best understood in terms of “Fluttering”

While the **Brain Dances** on the **Cusp** of **Shifts** the EEG “Flutter” just the same as the **Wings of Hummingbirds**

17



DYNAMIC SYSTEM CONTROL PROCEDURES

Extrinsic Constraints- Anything that limits the range of variability in nonlinear systems

Chaotic Perturbation- Randomly changing individual controls (like Linear systems do)

Entrain & Migrate -Sound systems- ‘Listening System’, Binaural Beats etc....

Synchronization through Chaos- The NeurOptimal® process

18

The Zengar **Basic Certification Guide**

The Zengar **Basic Certification Guide**




NEUROFEEDBACK MEASURES A SUMMATION OF ELECTRICAL CORTICAL ACTIVITY

NeuroOptimal® uses the 10/20 System and at C3 and C4 sites (Central Zone)

Black clip = Ground
Blue clip = Reference
Yellow flat Sensor = Active

19




NEUROFEEDBACK MEASURES - WAVEFORMS

We are NOT measuring EEG electrical activity PER frequency level, we are measuring a phasor of the constantly changing voltage of the EEG: the brain does not contain "waves" when measure on the scalp.

- * **Frequency:** "Windings" Per Unit Time
- * **Amplitude:** Radial Length Per Unit Time
- * **Phase:** Radial Rotation Per Unit Time

20



CARE MODEL

Comprehensive
Adaptive
Re-normalization
of the EEG

BY ZENGAR

21

The Zengar Basic Certification Guide



THE ROLE OF THE TRAINER...



- To be a **Detective**
- To **Empower** people to look within
- To help **Identify** any 'Shifts' or changes that have occurred
- To **provide** information
- To **Educate** people to the process and the **PASS** and **Community Resources** available to them

22




TRAINERS



- **Releasing Expectations** (not introducing expectations)
- **Educating Clients** (courses opportunities, using only our philosophy/language and branding)
- **Answering questions** (create answers ahead of time, so you are comfortable)
- Helping see changes (**monitoring** the shifts & DIFS)
- **Be ok if Client isn't "happy"**

23



TRACKING YOUR SHIFTS

Client Instructions:

Choose items that you circled on the "Tracking Your Shifts" The One Hundred, and mark them down on this sheet.

Then assign them DIFS.

Duration (how long that item lasts when it happens), **Intensity** (how intense it is when it happens, on a scale of 0-10) and **Frequency** (how often this item happens per day or per week) for each item.

In the following example, we will use feeling Frustrated as the item we wish to track.


Duration: When you feel frustrated, how long does it last? Is it a few minutes? A few hours? Estimate this, not counting when you are asleep and write it down for each item.

Intensity: What is the average intensity, 0-10 of your Frustration level when it happens?

Frequency: Using the past week as a guide, ask yourself "How often was I frustrated?" You can write down how many times you experienced something (3, for e.g.), but if you have a lot of something it is harder to ask "How many days this week did I feel angry?" So if you felt angry every day last week, your number would be 7/7/7 (and not 4/7). If something happens really often, measure how many times per day it is happening. It can also be interesting to log what triggers the item and how you deal with it when it happens.

Trainer Instructions:

On this sheet you will simply double check if the DIFS (Duration, Intensity, Frequency) have actually been assigned to each item the Client wrote down. If they have not, ask questions.



COURSES:


We offer a variety of different courses and certifications that go more into depth about the DIFS and other topics. Visit www.neurooptimal.com to learn more about our Basic, Advanced, Business and Re-Certification courses.

Tracking Shifts:

- Checklist 100
- WishList
- Pre- Session Evaluation

24

The Zengar Basic Certification Guide



NEUROOPTIMAL

Myth Busters

Rumor Has It...

Funny how rumors get started...One person thinks they understand how something works and says, "it must work like this" while another says "no, I heard it worked like this", and yet another says, "in my experience it must be like this". Before you know it suddenly, all these funny rumors, inventions and in some cases superstitions, have started and with little to no truth to it! So, I am going to take a moment to dispel a few rumors.

- 1 Allow at least ten nights sleep between NeuroOptimal sessions.
- 2 You must ease into the process by starting with the Initial session and eventually graduating to Regular as a process.

25



NEUROOPTIMAL

Myth Busters

Rumor Has It...

Funny how rumors get started...One person thinks they understand how something works and says, "it must work like this" while another says "no, I heard it worked like this", and yet another says, "in my experience it must be like this". Before you know it suddenly, all these funny rumors, inventions and in some cases superstitions, have started and with little to no truth to it! So, I am going to take a moment to dispel a few rumors.

- 1 Allow at least ten nights sleep between NeuroOptimal sessions.
- 2 You must ease into the process by starting with the Initial session and eventually graduating to Regular as a process.
- 3 If the Spectrograph has a predominant blue color, then add Omega 3 Oil or if the Spectrograph has 'red dots' within it, stop drinking caffeine.
- 4 Migraine headaches need 20 NeuroOptimal sessions, while ADHD needs 40 sessions.
- 5 I must have a QEEG to know where my brain is starting and what kind of progress I am making.
- 6 My Trainer must understand my diagnosis to best treat me with NeuroOptimal.

FALSE! NeuroOptimal is so safe that sessions can even be run twice per day if you want. Because it doesn't push the brain, NeuroOptimal cannot hurt you nor create side effects.

FALSE! Start with "Regular" and stay with "Regular", no matter who you are training. Since AutoNAV came to be, easing into your sessions using Initial, 2nd & 3rd is no longer necessary.

FALSE! The Spectrograph does not convey any information about outside influences such as these. All it does is give us a little information about the activity that occurred during the 30 second baseline. The colors mean nothing.

FALSE! Migraine headaches need 20 NeuroOptimal sessions, while ADHD needs 40 sessions.

FALSE! NeuroOptimal is a training, not a treatment. There are no special buttons we press because somebody has x, y or z. All it does is provide information to the brain about what it is doing, allowing it to make better decisions for itself. There is no way anyone can predict what shifts someone will experience, how big of a shift it will be or when that shift will happen.

FALSE! I must have a QEEG to know where my brain is starting and what kind of progress I am making.


FALSE! NeuroOptimal is a Dynamical Neurofeedback system that measures information while it is in movement, not in a linear moment in time. It is important that the measure is fluid and

FALSE! NeuroOptimal is a training, not a treatment. There are no special buttons we press because somebody has x, y or z. All it does is provide information to the brain about what it is doing, allowing it to make better decisions for itself. There is no way anyone can predict what shifts someone will experience, how big of a shift it will be or when that shift will happen.

FALSE! My Trainer must understand my diagnosis to best treat me with NeuroOptimal.

FALSE! NeuroOptimal is diagnostically agnostic. It is a mathematically driven program that does not need to know or understand anything about

26



NEUROOPTIMAL

Myth Busters

Rumor Has It...

Funny how rumors get started...One person thinks they understand how something works and says, "it must work like this" while another says "no, I heard it worked like this", and yet another says, "in my experience it must be like this". Before you know it suddenly, all these funny rumors, inventions and in some cases superstitions, have started and with little to no truth to it! So, I am going to take a moment to dispel a few rumors.

- 7 NeuroOptimal results (data) must be interpreted by an expert.
- 8 The NeuroOptimal Trainer must manipulate the software to create physiological changes in the CNS (Central Nervous System).
- 9 Zen 2 during a session can create more Dopamine, while Zen 3 will encourage more Serotonin production.
- 10 Zen 3 is too strong for this person or this age (for whatever reason).
- 11 I must stay with my client as
- 12 I believe my vibration bed.

FALSE! The information derived from NeuroOptimal is not necessary for continued training and is therefore not used in monitoring progress. Instead we use the 'progress tracking tools' to measure the DIFs (Duration, Intensity, Frequency and Shift) of the changes.

FALSE! The NeuroOptimal Trainer must manipulate the software to create physiological changes in the CNS (Central Nervous System).

FALSE! While each zen mode has a very precise 'mathematical' timing which results in each providing powerful feedback in its own way, there is no need to adjust the default timings. Since the introduction of autoNAV all the zen modes are fluidly navigated by the brain. Also remember, NeuroOptimal does not push the brain to 'do anything', no matter which mode you are in! Select "Regular" and trust NeuroOptimal to do its job.

FALSE! Zen 2 during a session can create more Dopamine, while Zen 3 will encourage more Serotonin production.

FALSE! Zen 3 is too strong for this person or this age (for whatever reason).

FALSE! I must stay with my client as

FALSE! I believe my vibration bed.

27

The Zengar Basic Certification Guide

NEUROPTIMAL
ADVANCED BRAIN TRAINING SYSTEM

FALSE! We do not need to 'help' clients through their process' by sitting with them. In fact, our being in the room may actually add an 'extrinsic constraint' to their process, which could slow it down. Let 'NeuroOptimal' and the individual dance together and get out of their way.

FALSE! We do not need to combine anything with NeuroOptimal® to make it effective! Sadly, although these other devices might be good products, they may slow down the NeuroOptimal® process due to their constraining effect on the brain.

13 Regular session is too strong/long for children. They must do shorter sessions until at least x years old.

FALSE! Do regular sessions all the time, no matter what age the person. If you find it a challenge to keep a child seated during a full session, quiet activities like reading, homework, or board games are fine. They can even watch a movie.

14 NeuroOptimal® Dynamical Neurofeedback™ system is no different than any other Neurofeedback approach.

FALSE! NeuroOptimal® is an approach to Neurofeedback that mathematically offers feedback by using a 4-D (4-dimensional) process of Duration, Intensity, Frequency &...

15 NeuroOptimal® creates neuroplasticity in the brain. NeuroOptimal® creates neurons in the brain. NeuroOptimal® re-connects neurons.

FALSE! The brain is designed to process information, and NeuroOptimal® offers exactly that – information. So there is no need to invoke physical effects of any kind. While people may have their...

28

Introducing **NeuroOptimal®** 3 The evolution of transformation

04.17.2018

NO3 DYNAMICAL NEUROFEEDBACK™

"In 1970 I had a vision of what brain training could be in the future: a simple, powerful, and safe way to relieve suffering that is so easy and intuitive that anyone can use it. NO3 is that elegant vision and it's here now"

- Valdeane W. Brown, Ph.D - Founder

29


Test length (seconds)
60

Start Test Stop Test Clear Chart Left Right

GENERAL TOUR NO3 VIDEO

NeuroOptimal® 3 Tutorials:
General Tour

30




31

NEUROPTIMAL
ADVANCED BRAIN TRAINING SYSTEMS
TRAINING YOUR MIND TO BE YOUR BEST


The Waterfalls: Introduction

- The **Waterfalls Display** replaces the Legacy **Spectrograms**
- Multiple **Perspective Views** accessed by a **Tap** with their own **Grab and Drag, Zoom and Reorient View**
- **Line Noise** is removed by default
- You can **Zoom by Intensity** and you can set different **Regions of Interest (ROI)** for greater precision and clarity
- You can migrate **ROIs** between the Left and Right Channel **Waterfalls**

32




33




Cross Ambiguity Function

- The **Cross Ambiguity Function (CAF)** replaces the **Legacy CCACs**
- CAF** is the larger group of mathematical algorithms. The **CCAC** is a particular way of determining the **CAF** for the transition of **Eyes Open to Eyes Closed**
- NO3** calculates **CAFs** from 3 minute segments of **EEG** throughout a **Training Session**
- These segments overlap 1 minute, making the **CAFs** even more robust

34



After Session Views



35



The Path of Self-Empowerment

• Website, PASS Membership & Resources

HELP Solve issues at any step along the way

Before a problem arises: Review the help guides, manuals, and resources. If you still have a problem, contact the Technical Support team.


1. Visit the website and review the help guides, manuals, and resources.
2. Contact the Technical Support team via email or phone.
3. Visit the website and review the help guides, manuals, and resources.
4. Contact the Technical Support team via email or phone.
5. Visit the website and review the help guides, manuals, and resources.
6. Contact the Technical Support team via email or phone.
7. Visit the website and review the help guides, manuals, and resources.
8. Contact the Technical Support team via email or phone.
9. Visit the website and review the help guides, manuals, and resources.
10. Contact the Technical Support team via email or phone.

How to Help-Yourself...

The Guides, Manuals, PDF's and other Resources for Self-Empowerment

36


The Zengar **Basic Certification Guide**



37



38




39









The Zengar Basic Certification Guide



43




44




45


The Zengar Basic Certification Guide




46



47



48



49

ZENGAR
INSTITUTE®

The Zengar **Basic Certification Guide**
