

SYSTEM MAINTENANCE

Zengar[®] recommends performing the following maintenance on your system routinely and in preparation for a software update or ZenConnect:

- 1. Set aside enough time to complete the process**
Plan to remain with the system during the process to allow any prompts that occur
- 2. Confirm at least 20GB is free and available on your SSD/HD**
Custom media files can be moved to external drive/microSD card if available
https://neuroptimal.com/technical-support-faqs/#check_my_storage
- 3. Check for and complete all Windows 10 updates**
<https://neuroptimal.com/technical-support-faqs/#faq4>
Do not upgrade to Windows 11 at this time
- 4. Run Disk Cleanup**
<https://neuroptimal.com/technical-support-faqs/#faq5>
- 5. Optimize the SSD/HD**
<https://neuroptimal.com/technical-support-faqs/#faq6>
- 6. Reboot the system**
- 7. Initiate the ZenConnect or update process if desired**
ZenConnect is not required for the update process
Never stop an update in progress