




How to find the System ID on your NeurOptimal® system.

Zengar Institute Inc.


1. In the Windows Cortana Search field, type **About your PC**.



2. Click the **About your PC** icon.



3. In the About Window, under **Device Specifications** you will find the Device name for your NeurOptimal® system.



About

Your PC is monitored and protected.

- ✓ Virus & Threat Protection
- ✓ Firewall & Network Protection
- ✓ App & browser control
- ✓ Account protection
- ✓ Device security

[See details in Windows Security](#)

Device specifications

Device name	JSmith_01
Processor	Intel(R) Core(TM) i5-7300U CPU @ 2.60GHz 2.71 G
Installed RAM	8.00 GB
Device ID	ACD92B13-BEC6-4B0F-A430-DAE8939FD926
Product ID	00330-62945-24897-AAOEM
System type	64-bit operating system, x64-based processor
Pen and touch	Pen and touch support with 10 touch points

Rename this PC

Windows specifications

Edition	Windows 10 Pro
Version	19H2